

HARD BODY ARMOUR PLATES


Data sheet TenCate Multi-Light™ CS-751 SA (Torso)

ADVANCED ARMOUR

CS-751 SA (TORSO)









| | |
|---------------------------------|---|
| Configuration | Stand alone |
| Type of protective cover | Black textile or TenCate Rosist® |
| Product application | Tactical Body Armour, External and Uniform Carriers |

SPECIAL THREATS TESTED

| | |
|--|---|
| VPAM 9 |  |
| 7.62 x 39 API BZ | |
| 5.56 x 45 SS109 | |
| 7.62 x 39 AK 47 MSC | |
| 7.62 x 51 M80 | |
| <p>The above list of ammunition is not exhaustive and may not include threats that can be defeated by this product</p> | |



SPECIFICATIONS

| | | | |
|---|---|--|---|
|  |  |  |  |
| Rifle ammo resistant | Multi-hit capability | Stop in plate capability | Lightweight |
| 7.62 x 39 API BZ | 6 shots | Against listed threats | 3.25 kg ± 0.1 kg for textile cover 3.3 kg ± 0.1 kg for Rosist® |
|  |  |  |  |
| Dimensions | Shape | Thickness | Composition |
| 250 x 300 mm ± 2 mm | Multi-curved | 19 mm ± 2 mm for textile cover 20 mm ± 2 mm for Rosist® | Monolithic Ceramic HPPE backing |

HANDLING

- › Do not sit on the insert
- › Do not drop the insert on the floor and avoid other hard shock(s)
- › Do a visual inspection of the insert every six-month
- › In case of visual damages, send the insert back for inspection / repair
- › Do only clean with pure water or mild soap

GUARANTEE PERIOD

The plate is expected to preserve its conformity with the specifications for a period of at least 5 years during storage from the date of delivery, provided that:

- › The product is stored under dry conditions at normal ambient temperature
- › The product is stored vertically
- › The product has not been subjected to peak temperature exceeding 70°C
- › The product is protected against long term direct exposure to sun light